

Multi-function 7 in 1  
Skin Analyzer Facial  
smart ice blue skin  
management system

For skin analysis skin cleaning and skin care










For skin analysis skin  
cleaning and skin care

# Product details



-  Skin analysis pen
-  Ion bubble pen
-  High frequency skin scrubber
-  Gold RF pen
-  Ice head



-  Spray gun pen
-  Ultrasonic pen
-  Touch screen
-  ABC Solution bottle

Technique Creates Beauty

A  
S  
I  
S  
E  
R  
V  
I  
C  
E  
A  
N  
D  
B  
E  
S  
T  
V  
A  
L  
U  
E

## Handles introductions



### Skin analysis pen

Panoramic optical intelligent skin detection-high-definition imaging technology Combined with dermatologist evaluation + medical big data sampling + AI intelligent diagnosis and analysis, the accuracy rate is up to 98%. Subdivide 8 categories of skin conditions to detect skin microscopically.

### Ion bubble pen

Air pump imported from South Korea, with strong suction power, can suck out deep pore dirt, combined with a panoramic skin macro detection system, clean pore dirt without dead angle.



### Spray gun pen

Nano-level ultra-micro atomization, decomposes the essence and skin care ingredients into nano-level small molecules, and penetrates deep into the muscle bottom through atomization and penetration technology to provide deep nutrition

### Ultrasonic pen

1 million high-frequency ultrasonic waves per second, reaching the depth of the fascia layer. Turn skin care products into diffused into deep layers, reaching the bottom of the muscle. Deep diffusion skin care product absorption



### Ice head

Semiconductor cooling + air cooling + water cooling technology can reduce the temperature to -5 degrees, feeling comfortable, soothes and repairs sensitive skin, shrinks pores and locks moisture and nutrients

### Gold RF pen

The 24K gold-plated radio frequency head fits the skin to help skin care products absorb and focus on heat directly to the bottom of the muscle, stimulate the proliferation of collagen and rebuild the collagen fiber network.





### High frequency skin scrubber

Ultrasonic 24KHZ-25KHZ high-frequency vibration, impacts the dirt in the deep pores, will blackheads, small dirt such as acne and make-up residue will be removed from the skin surface without damaging the skin barrier.

## Function of machine



1. Facial lift: skin rejuvenation, fat removal
2. Eyes beauty: eyes pouch removal, eyes dark circle removal, wrinkle around eyes removal
4. Whitening and anti-wrinkles
5. Take away the molecular structure of cells; skin completely open power operation in the skin surface infiltration head with a method of grain is slightly along the skin and smooth movement of extrusion.
6. Clean the skin, daub hydrating moist products
7. Cool frozen, effect on skin, lock moisture and nutrients.



Technique Creates Beauty

# Software Design



Skin treatment menu



Homepage



- Immediate visible results
- Pleasant and soothing
- No downtime
- Safe for all skin types
- Optimal cleansing and exfoliation
- Maximum infusion of active ingredients
- Superior active nourishment

Technique Creates Beauty

## ■ Technique Specification

Vacuum range	700 mmHg
Output control	Touch screen & handle switch
Gold RF	Rapid response multipole (energy level 10 regulation)
Electrophoresis	1MHZ( energy 1 level 0)
Ultrasonic ultrasound	1MHZ( energy 1 level 0)
Frozen head	Water cooling & air cooling (energy level 10 regulation)
High-frequency scalpel	25KHZ( energy level 10 regulation)
Ion Shotgun	mmHg 2000
Interface	10.1 inch computer touch screen
Input voltage	110-240 V 50-60HZ 150AV