





MEDI-TECH SOLUTIONS (The Complete Aesthetics)



DIDDE-LASER

Technical Specifications:

- Wavelength Laser Power Pulse Width Frequency Energy Fluence Skin Cooling LCD Screen Input VAC Dimension Net Weight
- Triple Wave (755/808/1064/nm) 1000W 5~160ms 1~10Hz 0~100J/cm2 2TEC+Sapphire 8.4 inch Touch LCD 110/220VAC / 50-60Hz 520x550x1250cm 15kg

DIDDE LASER FOR PERMANENT HAIR REDUCTION

Feature:

- Long Life time upto 50 million shots
- Cooling Temperature –10°c to 5°c
- No degradation of Energy.
- Bigger Spot Size 14x23mm
- Pain less Technology



Q-SWITCHED ND YAG LASER

Technical Specifications :

Laser type	Q-switch Nd: YAG Laser
Wavelength	1064nm, 532nm and 1320nm
Pulse Energy	100-1000mj
Pulse Width	3ns
Facular Diameter	1-8mm
Pulse Frequency	1-10 Hz (adjustable)
Screen	8.4" True Color Touch Screen
Electrical Requirements	AC100-240V, 20A max., 50/60Hz
Dimension	38*32*48cm (W*D*H)
Net Weight	25 kg

Application

- Tatto Removal
- Freckle Removal, Nevus Ota Removal
- Birthmark Removal, Embroider eyebrow Removal
- Skin rejuventation, Whiten, acne, tighten pores (Carbon tip needed)



D

MTS

VERTICAL Q-SWITCHED ND YAG LASER

Technical Specifications

Laser type Wavelength LCD screen Single-pulse Energy Double pulse energy Power Pulse width

Frequency Spot size Dimension of machine N.W. Voltage E-optical Q-switched Nd YAG laser 1064nm&532nm&1064nm static light 10.4' TFT 1.3J (1064nm) ; 650mJ (532nm) 2.6J (1064nm) ; 1300mJ (532nm) 2000W 6ns @ 1064nm/532nm Q-switched 300us@1064nm Quasi-long pulse 1-10 2-10mm adjustable 92x34x114.5 cm 98kg 220V / 110V

Application :

- Pigment deposit dispelling, even effective for chlosma/melasma
- Color tattoo removal
- Eliminates the age pigment, spot, birthmark
- Pigment changes
- Skin rejuvenation, skin toning
- Lip line and eyes line removal

Renuva Non Ablative Fractional Laser System

Technical Specifications

Wavelength Pulse Energy Beam density Tips Mode Weight Dimensions 1550nm Up to 70mj Up to 500 micro-beams/cm² 5 Fractional Scanning 20kg 492mmx400mmx353mm

Features :

- Skin Rejuvenation (photodamage, rhytides, skin texture.)
- Wrinkle Removal (such as aging wrinkles, large pores, rough skin, etc)
- Scar Repair (such as acne, surgical scar, caesarean scar, burn scar, etc)
- Vergeture Repair
- Reduce Pigment (such as dark circles etc)
- Hair Regrowth



HAIR REDUCTION & SKIN CARE

Technical Specifications:

IPL SHR ELIGHT

Power	2500W
Capacitors	3*12000uf
Energy Density	1-25j/cm2 (Europe safety standard)
Frequency	1-10 shots/seond
Wavelength	650nm-950nm & 570nm-950nm
Hand piece	2; SHR 15*50mm2 & SSR 12*33 mm
Pulse Duration	1-15ms
Cooling system	Sapphire+air+water+SPT cooling
Electrical Requirements	100-240VAC,20A max, 50/60Hz

Application :

- Hair Reduction
- Skin Rejuvenation, Pigmentation, Acne, Vascular Therapy, Wrinkle Removal, Face Lifting

CO2 FRACTIONAL LASER

Technical Specifications

40W RF Metal Tube
10600nm
0.2mm
Fractional, Ultrapulse and Vaginal
7, Circle /Triangle/Square/
Rectangle Rhombus/Line/Ellipse
7 joint arm
Max 20mmx20mm Air cooling 110V/220V, 60/50Hz

HIFU SKIN TIGHTENING SYSTEM

Technical Specifications

Туре
Frequency
Power
Energy
shots Gap
Power Source

High Intensity Focus Ultrasound 4MHz 200W 0.1---2.0j 1.0---3mm AC 110V~220~50Hz/60Hz

Features :

- Patented SPT, FCA technology, focus on dark and light hair.
- Germany imported xenon lamp, stable, effective
- 10 shots/seond fast hair removalThermage (Fractional RF)
- handpiece optional
- -2°C Speedy cooler say goodbye to pain





Application :

- Wrinkle removal & skin resurfacing
- Pigmentation (Freckle, Chloasma, Sun damage)
- Acne removal, Acne Scar
- Scar removal, Stretch Mark
- Cutting Surgical (Warts, Nevus etc)
- Vaginal tighten, Vulva rejuvenation (VRL function optional)



Application :

HIFU is Non Surgical Face Lift Treatment - HIFU (High Intensity Focused Ultrasound) Technology, the aging effects will show on the skin like the wrinkles and the eye-bags. It is the signals to alert you your collagen and elastin are losing from your body. The collagen is the sky structural component to support your skin.



Cryolipolysis

Technical Specifications		
Power	1500W	
Focused Technology	360 CryoSculpt+Vacuum	
Handpiece	PLUS : 145x135x50mm HPI : 200x50x70mm HPII : 150x45x50mm HPIII : 120x45x50mm MINI : 65x20x20mm	
Cooling Stand-by Working	-10°C~0°C Continuous contact cooling	
Vacuum pressure	Level1~10 adjustable	
Display	8" True Color LCD Touch Screen	
Net Weight	36 KG	
Dimensions	575x570x1200mm	

(E 🎱 🛞 🛞 FDA

Application :

- Face thinner, skin tightening
- Fat reduction, cellulite reduction
- Body slimming, body shaping
- To promote and accelerate the body's metabolism

Features

- Treat 2 Areas Simultaneously in only 40 Minutes.
- Upgraded Technology for Enhanced Results.

MEDICAL CE APPROVED





MICRO NEEDLE RF & FRACTIONAL RF

Technical Specifications

System Type	Mono & Bipolar RF/ Bipolar RF
Frequency	2MHz
Total RF Treatment Level	Up to Level 10
Dual Operation Mode	FX RF, MN FX RF Sub-lative
Max. Energy	50W
Depth of Needle	0.5-3mm (Adjustable)
Display Switch	10.4 Inch LCD Foot Switch
Dimension (mm)	280 mm(w)x505mm (L)x1100mm(H)
Electrical Requirements	AC/100-240 V, 50-60 Hz

Features :

Pinxel is designed specially for delivering RF energy directly to epidermis or dermis layer in order to stimulate and remodel collagen effectively by either RF micro needing mode (invasive) or RF matrix mode (non-invasive). It can be used for both epidermal applications and dermal application in one platform without changing any operational environment.



SYNERGY M340

Technical Specifications

Model No.
Radio Frequency Powe
Focused Technology
Delivery Type
RF Frequency
Pump Suction Rate
Infrared Lamp
Light Spectrum

Handpiece

Cooling Display Net Weight Dimensions MED-340 Up to 50W Vacuum+RF+Roller+LED+IR Pulse;CW 10MHz 25cub.m/h up to 20W 640nm (infrared); 690nm (LED) 1xVacuum+RF (56x80.5mm) 1xVacuum+RF(34.5x48mm) 1xVacuum+RF(18x25mm) 1xVacuum+RF(18x15mm) Self-contained 10.4" True Color LCD Touch Screen 51kg 56x47x110cm

Application :

- Versatile 5 technologies in 1: High Efficiency
- Handle self-owned operation screen: Convenient
- Ergonomic design : Easy to handle
- Convenient Operation with handle button 4 different rolling directions with multi-dimensional
- effects High precision processed treatment head
- Non-invasive delicate treatment for Face Carving
- Uniform RF energy output brings accurate treatment
 - results
- Nonperishable Comfortable Medical Materials Head



Permanent Tatooing / micro pigment unit

Feature:

- 1. Permanent makeup application:
- Lip tattooing(out line and full lip)
- Eveliner tattooing

Sneed.

Motor:

Color:

- Breast areola beauty
- Needle length: 28000 cycles/min Speed of motor: 10 levels of speed Imported for Swiss white
- 2. Micro-needling application:
- Skin tightening, Lifting, Rejuvenation
- Eyebrow tattooing (Can make very thin) Cure acne scar & wound healing
 - Improve wrinkles and fine lines



Advantage:

- Two handpiece can work at same time, speed can be adjusted independently
- Standard motor for micropigmentation, powerful and strong
- Specially tailored needles for PMU, steady movement with no shaking

PRO MEDI FACIAL

Technical Specific

Rated Power	110W
Micro - current power	10W
Vaccum pressure	0-75Кра
Machine Material	ABS
Handpiece Number	1. Liquiddermabrasion handle with 8 tips
	 Diamond dermabrasion handlic with 9 Tips
	3. Multipolar RF head
Dimension of Machine Net Weight	35x29x40cm 7.5kg
Voltage	AC220V±20% ,50Hz AC110V±20% ,60Hz
	ACTIOVIZO/0,00HZ

Indications

- 1. Deep cleaning of the skin Clean blackheads, acne, mites
- 2. Shrink pores Improve the dark skin brighten skin tone
- 3. Moisturizing Improve skin quality
- 4. Photoelectric device Preoperative and postoperative care



CE 🖉 🛞 🛞 🗗 Medical CE Approved



DIGITAL MESO GUN

Features :

- Non-invasive, No blood, No Pain, Effective with no downtime! •
- Powerful electroporation energy!
- Multifunctional, for Skincare and Hair Care with two different cartridges!
- Allows a wide variety of active ingredients to be introduced into the skin!
- Ergonomically designed applicator, Super lights!
- Match with match various types of syringe!





Parameter

Power Supply Adapter Vibrating Frequency can use 8 hours after charging speed . Needle length Needle No.:

Used Method

Different needle length takes a different using method The Size of 0--0.50mm use method is: in the skin on the circle treatment. The size of 1.00--2.00mm use method is : likes the stamp roller using way stamp on the skin)

Rechargeable 4.2v-500MA 3000, 4500, 7500r/m

3 level vibration speed controlled from 0.1mm-2.0mm adjustable 9/12/21 pin needle, nano needle tips

Indication :

- Ilication : Winkle removal, improve scar, skin rejuvenation For the narrow focus of the site can be treated (glabella, depth, winkle, scar) Shrink pores & Lifting Micro pigmentation & cosmetic tattooing

ELECTRO PORESION NEEDLING (EPN)

- The part of the Needle Set to be attached You can change the multi-needle according to types of skin and face.
 - LCD Display Easy menu control with indicator light.

Β.

- Speed
 : Speed of the motor, Adjustable from level 1 to Level 3

 EP:
 EP Level, Adjustable from Level 0 to Level

 Time:
 Duration of use, up-counting 30Min
- Control Panel Menu control button / Adjustable operation mode / Adjustable Level / Adjustable Time Reset
- Start, Stop You can adjust only the level of speed during opera You cannot adjust MODE and TIMER during opera
- The adjusting part of the needle penetration depth It enables you to offer customized treatments by adjusting the needle penetration depth.

MDAUNIT

MICRO DERMABRASION SYSTEM

Features :

- Desktop Model
- Crystal /Diamond Mode
- Uniform Abrasion
- Maintenance Free

MDA PEN

Features :

- Portable Model
- Diamond Mode
- Digital Display
- Handy to Use



The therapy exfoliates the top most layer of the skin (stratum conium) It creates a "closed system" creating a vacuum which suck in the dead skin and debris. Microdermabrasion therapies are used to reduce scarring, lightening hyper pigmentation and making the skin lighter, younger, and smoother.

CAUTRY PEN WITH ELECTRO PORESION TIP

Condition Material: Plastic, Color : White Package Weight Item Length Function

: 100% Brand New Type:

- : Approx. 261g ,
- : Approx.17.6cm / 6.9inch
- : Can be applied to spots,
- dot mole, tattoo removal

Features :

- 1. Adopt the latest technology, it is safe, comfortable and fast to remove the freckle and moles without any side effect.
- 2. Use 3 adjustable scan spot power output for different treatment, which is reliable, useful and convenient.
- 3. Unique electric ion technology makes it remove spots without hurting the normal skin, no bleeding and little pain.

CAUTRY PEN



.New launch Product Image....

Technical Specifications

Xecl Lamp
308nm
30mW/cm ² ±20%
30-4200mJ/cm ²
Max. 16cm²
1s~140s
3s/6s/9s/12s/15s/18s/21s/24s/27s
Continuous
27x24x26cm
100-240V, 47-63Hz

Feature :

- The 308 nm single wavelength excimer light output via the medium XeCl
- High dose, highly effective and fast result
- Targeted treatments, maximum protection of healthy skin
- Up to 16cm² spot size for short treatment times of larger areas, fast treatment
- Use-friendly Interface, convenient and easy to operate









308NM EXCIMER SYSTEM



The Protable, Highly Effective Solution for Monochromatic UVB Therapy Psoriasis Vitiligo Eczema Dermatitis



for more Informations:

() : +91-93112 87626, 011-4704 3238

- meditechsolutions.mts@gmail.com
 info@meditechsolutions.in
- 🕀 : www.meditechsolutions.in
- S: 7827263140