

# Q SWITCHED Nd YAG LASER TATTOO REMOVAL



**MEDICAL CE APPROVED**



## How it Works

If you're worried about lengthy sessions with a laser, don't be. Q switched laser with very short pulse durations 3ns, which help keep skin free of scarring. To further protect the skin, treatments are usually done over several sessions. The number of sessions you need will depend on your tattoo. Every tattoo is different and they come in all shapes, sizes and colors, meaning some tattoos are easier to remove than others. Professional tattoos (one you get at a tattoo parlour) tend to take longer to remove than those done by amateurs (one you get at a party after one too many tequila shots).





## Applications

- Tattoo Removal;
- Speckle Removal; Nevus of Ota Removal;
- Birthmark Removal; Embroider eyebrow Removal;
- Skin rejuvenation, whiten, acne, tighten pores  
(Carbon tip needed)

## Features

1. Up to 10Hz for faster tattoo removal
2. 8" Color touch screen easy operation
3. Infrared aiming beam for accurate operating
4. Lamp life time over 1,000,000 shots
5. ABS housing for quick maintenance work



**1064nm** head for black, blue, brown tattoo ;  
**532nm** head for red, orange, pink and colorful tattoo.



**1320nm** Carbon peeling head optional, for Skin rejuvenation, whiten, tighten pores, acne

