

MEDI-TECH SOLUTIONS

Q Switch ND YAG Laser Machine

m

b

Z

m







Q Switch ND YAG Tattoo & Pigmentation Removal Machine

Application



MEDI-TECH SOLUTIONS

Advantages

01		25	0
RC Carlot	-		
		·	
	Statutes.		-
for Australiants	and hear to pe		
*****	 Bernards, Fill, and Roberts, Roberts, March Science, Roberts, W. Science, Science, W. Science, Science, W. 		2
1.1		and the second	-
1000 C			
-	18.	7	

1. TUV Medical CE approved

Strict quality testing process, safe treatment and stable working.

Þ

C

00

P

C

IT.



2. Aiming light accurate operating Aiming light can be turn on/off, longer lifetime



3. High-efficiency Dermal Pigment & Tattoo Removal Double energy handpiece especially focus on dermal pigment and colorful tattoo removal, permanent results and less sessions is needed.



4. Super long lifetime Lifetime of xenon lamp more than 1,000,000 shots.



5. Intelligent interface Five different tips for various treatments can be selected

MEDI-TECH SOLUTIONS

Technique Specification

wavelength	1064nm, 532nm, 1320nm, Fractional,		
	Focus Array; Five Tips Standard		
Laser Power	500W		
Energy	300-2000mj		
Yag bar size	Φ7+φ8		
Display	8.4" True Color LCD Touch Screen		
Cooling	Air+Water+Twin Radiators Cooling System		
Net Weight	25KG		
Dimensions (WxLxH)	40*36*46cm		



Company

Name: MEDI-TECH SOLUTIONS Email : meditechsolutions.mts@gmail.com